
Type	Component	Issue ID	Description
Improvement	Visit Management	ASTHLS-645	The visit dialog window now has a scroll bar if the number of available visits/orders exceed the number of available rows.
Improvement	Visit Management	ASTHLS-647	The (pre)-selection of the visit/order in the visit dialog can now be configured via defined criterias. The default selection remains the last received one.
Improvement	Visit Management	ASTHLS-648	It is now possible to activate via a property a warning if a user select a visit/order to be performed in future. The number of days in future can be configured (default is one day)
Improvement	Visit Management	ASTHLS-649	A new field to store visit type has been added. This value can be optionally seen in the visit dialog (activated by property)